



Chem 1112 Spectrometry Dye Solution Preparation	Identifier: CSP-0017
	Revision: 0
	Page: 1 of 4

ISU Chemistry Department	Stockroom Procedure	Effective Date: 10/01/2019
--------------------------	---------------------	----------------------------

1. INTRODUCTION

This procedure provides guidance in the preparation of dyed water solutions for the CHEM 1112 spectrometry laboratory experiment.

2. PRECAUTIONS AND LIMITATIONS

- 2.1. Food coloring is able to stain clothes
- 2.2. Orient the cuvette with its frosted sides facing the front and rear of the instrument.

3. APPARATUS AND MATERIALS

- 3.1. Spectronic™ 200 Spectrophotometer
- 3.2. Plastic cuvettes
- 3.3. 2L beaker
- 3.4. 4L beaker
- 3.5. 4L carboy
- 3.6. Glass stir rod
- 3.7. Transfer pipettes, any size
- 3.8. Waste beaker, around 150 mL
- 3.9. Wash bottles
- 3.10. Paper towels
- 3.11. Label maker
- 3.12. Clear packing tape

4. REAGENTS

- 4.1. Water
- 4.2. Food coloring

5. INSTRUCTIONS

- 5.1. Getting water
 - 5.1.1. Measure 2L of DI water in the 2L beaker.
 - 5.1.1.1. Pour the water into the 4L beaker.
 - 5.1.2. Repeat steps 5.1.1. and 5.1.1.
- 5.2. Dying the water



Chem 1112 Spectrometry Dye Solution Preparation	Identifier: CSP-0017 Revision: 0 Page: 2 of 4
--	--

ISU Chemistry Department	Stockroom Procedure	Effective Date: 10/01/2019
--------------------------	---------------------	----------------------------

- 5.2.1. Add approximately 50 drops of dye to the 4L of water.
 - 5.2.1.1. Record the exact amount of drops added.
- 5.2.2. Stir the water with a glass stirring rod until it is a homogenous solution.

5.3. Preparing the Spectrophotometer for operation

- 5.3.1. Turn on the spectrophotometer.
- 5.3.2. Remove any cuvette in the sample compartment.
- 5.3.3. Close the lid.
- 5.3.4. Wait for the instrument to finish its series of self-tests.

5.4. Select the correct settings

- 5.4.1. Select SPEC 200E Modern Interface.
- 5.4.2. Use the UP, DOWN, LEFT, RIGHT arrows to toggle between the options.
 - 5.4.2.1. Switch the Measurement Mode to Abs.
 - 5.4.2.2. Change Measurement λ to the proper wavelength with the knob.

Note: Every color absorbs at a different wavelength. In this procedure, λ for the blue dye is 631 nm, and λ for the yellow dye is 427 nm.

- 5.4.3. Select GO.
 - 5.4.3.1. Press Enter.

5.5. Running Samples

- 5.5.1. Blanking the instrument
 - 5.5.1.1. Place a cuvette with DI water in the sample stage.
 - 5.5.1.2. Press the AUTOZERO button.
 - 5.5.1.3. Wait for the autozero to finish.
- 5.5.2. Press ENTER to freeze data collection.
- 5.5.3. Open the sample compartment.
 - 5.5.3.1. Remove the blank cuvette from the sample chamber.
 - 5.5.3.2. Put it to the side.
- 5.5.4. With a transfer pipette, fill an empty cuvette with the dyed water.
 - 5.5.4.1. Put it in the sample compartment.
 - 5.5.4.2. Close the lid.



Chem 1112 Spectrometry Dye Solution Preparation	Identifier: CSP-0017 Revision: 0 Page: 3 of 4
--	--

ISU Chemistry Department	Stockroom Procedure	Effective Date: 10/01/2019
--------------------------	---------------------	----------------------------

- 5.5.5. Press ENTER to resume data collection.
- 5.5.6. Record the Abs next to the amount of food coloring added.

Note: *The ABS needs to be between 1.5 and 2 for both colors of water.*

5.6. Increasing the Abs value

- 5.6.1. Add drops of food coloring to the water.
 - 5.6.1.1. Record the new total drops added.
- 5.6.2. Repeat steps 5.5.4. - 5.5.6. until the desired abs is reached.

5.7. Decreasing the Abs value

- 5.7.1. Pour out 100mL of dyed water.
- 5.7.2. Add 100mL of pure DI water.
- 5.7.3. Repeat steps 5.5.4. - 5.5.6.
- 5.7.4. If the Abs still needs to be decreased, repeat steps 5.7.1. - 5.7.2.

5.8. Finishing

- 5.8.1. Transfer the dyed water with the correct Abs to a 4L carboy.
- 5.8.2. Label the carboy.
 - 5.8.2.1. Apply a piece of clear tape to the carboy.
 - 5.8.2.2. Print a label with the correct concentration (see 6.1.).
 - 5.8.2.3. Apply the label above the clear tape.
 - 5.8.2.4. Put an additional piece of tape over the label.

5.9. Cleanup

- 5.9.1. Rinsing cuvettes
 - 5.9.1.1. Squirt DI water into the cuvettes.
 - 5.9.1.2. Empty the cuvette into a waste beaker.
 - 5.9.1.3. Repeat 5.8.1.1. - 5.8.1.2. two times.
 - 5.9.1.4. Put the cuvette upside down on a paper towel.
- 5.9.2. Pour the contents of the waste beaker down the drain.
- 5.9.3. Rinse all glassware DI water three times
- 5.9.4. Throw all other used materials in the trash.



Chem 1112 Spectrometry Dye Solution Preparation	Identifier: CSP-0017 Revision: 0 Page: 4 of 4
--	--

ISU Chemistry Department	Stockroom Procedure	Effective Date: 10/01/2019
--------------------------	---------------------	----------------------------

6. Calculations

6.1. Beer's Law

$$Abs = \epsilon b c$$

Where:

ϵ = Molar absorptivity

b = path length = 1 cm

c = concentration in Molarity

Note: *The molar absorptivity constant changes for every solution. In this procedure, it is equal to 1.3×10^5 for the blue dye, and 2.73×10^4 for the yellow dye.*

$$c = \frac{Abs}{\epsilon b}$$